

---

## BCOR (Clone C10)

### Mouse Monoclonal Antibody

**Specificity:** Human. Others unknown

**Immunogen:** A synthetic peptide against amino acid 1-300 mapping at the N-terminus of BCOR of human origin


**Ig Class:** IgG<sub>1</sub>/κ

**Storage:** Store vial at 4°C. When stored at 2-8°C, this antibody is stable for 24 months

**Staining procedures:** Use formalin-fixed and paraffin-embedded sections. *Retrieval conditions:* Pretreatment of deparaffinized tissue with heat-induced epitope retrieval is recommended. *Detection methods:* Polymer anti-mouse/rabbit IgG detection system. *Working dilution:* 1:100 *Positive Control:* Endometrial stromal sarcoma. *Cellular Localization:* Nucleus. *Intended Use:* In vitro diagnosis (IVD).

**Description:** BCL-6 interacting corepressor (BCOR) is a 1,755 amino acid protein that associates with histone deacetylases to transcriptionally repress BCL-6. BCOR mRNA is ubiquitously expressed in various tissues. However, the BCOR protein expression in adult tissue is limited to testis only. Recently, the BCOR rearrangement was detected in a subset of small blue round, spindle cell, and undifferentiated sarcomas, which results in overexpression of BCOR protein in these tumors, including small round cell sarcoma (~100%), clear cell sarcoma of kidney (100%), high-grade endometrial stromal sarcoma (100%) with *YWHAE-NUTM2* translocation, and subset of synovial sarcoma (50%). The list of BCOR positive tumors is still expanding.

**Supplied As:** Purified antibody with 0.2% BSA and 15mM sodium azide.



*Formalin-fixed, paraffin-embedded human small blue round cell sarcoma stained with anti-BCOR antibody using peroxidase-conjugate and DAB chromogen. Note the nuclear staining of tumor cells*


**Cat. #Z2334 (1.0 ml)**

---

**ZETA** Corporation  
65 N 1<sup>st</sup> Ave, Ste 202C  
Arcadia, CA 91006, USA  
Tel: (626) 355-2053  
<http://www.zeta-corp.com>



Zeta Corporation



**Emergo Europe B.V.**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
Tel: +31 70 345 8570