

p63 (Clone ZM70)

Mouse Monoclonal Antibody

Specificity: Human, mouse, and rat. Others-not known.

Immunogen: Recombinant full-length human p63 protein


Ig Class: IgG2b/ κ

Storage: Store at 2-8°C for up to 2 years for concentrate form and 1 year for predilute form

Staining procedures: Use formalin-fixed and paraffin-embedded sections. *Retrieval conditions:* Pretreatment of deparaffinized tissue with heat-induced epitope retrieval is recommended. *Detection methods:* Polymer anti-mouse/rabbit Ig detection system. *Working dilution:* 1:100-200. *Positive control:* Skin, prostate, *Localization:* Nuclear. *Intend use:* In vitro diagnosis (IVD).

Description: p63 is a homolog of the tumor suppressor p53. It is identified in basal cells in the epithelial layers of a variety of tissues, including epidermis, cervix, urothelium, breast and prostate. p63 was detected in nuclei of the basal epithelium in normal prostate glands; however, it was not expressed in malignant tumors of the prostate. As a result, p63 has been reported as a useful marker for differentiating benign from malignant lesions in the prostate, particularly when used in combination with markers of high molecular weight cytokeratins and the prostate-specific marker AMACR (P504S). p63 has also been shown to be a sensitive marker for lung squamous cell carcinomas (SqCC), with a sensitivity of ~90%. Specificity for lung SqCC, vs. lung adenocarcinoma (LADC), is approximately 80%. In breast tissue, p63 has been identified in myoepithelial cells of normal ducts.

Supplied As: Purified antibody in PBS buffer with 0.2% BSA and 15mM sodium azide.



Formalin-Fixed, paraffin-embedded human prostate carcinoma stained with anti-p63 using peroxidase-conjugate and DAB chromogen. Note nuclear staining of basal cells in benign glands whereas carcinoma glands are negative

Cat. # Z2380

ZETA Corporation
65 N 1st Ave, Ste 202C
Arcadia, CA 91006, USA
Tel: (626) 355-2053
<http://www.zeta-corp.com>



Zeta Corporation

Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Tel: +31 70 345 8570