

PSA (ZR232)

Rabbit monoclonal Antibody

Host: Rabbit

Specificity: Human

Immunogen: Recombinant fragment (around aa150-250) of human KLK3 protein (exact sequence is proprietary)

Ig Class: IgG

Storage: Store at 2-8 °C for up to 2 years for the concentrate form and 18 months for the predilute form.

Specifications: PSA (ZR232) Rabbit monoclonal primary antibody detects PSA protein in formalin-fixed, paraffin-embedded samples by immunohistochemical (IHC) staining.

Staining procedures: Use formalin-fixed and paraffin-embedded sections. Retrieval conditions: Pretreatment of deparaffinized samples via heat-induced epitope retrieval is recommended. Detection method: Polymer anti-mouse/rabbit Ig detection system. Working dilution: 1:100-200. Positive Control: a sample known to contain PSA protein.

Localization: Cytoplasmic.

Intended Use: This antibody is intended for Analyte Specific Reagent (ASR) use only. PSA (ZR232) Rabbit monoclonal primary antibody is intended for laboratory use in the detection of adenovirus in formalin-fixed, paraffin-embedded tissue by immunohistochemical (IHC) staining. Qualified pathologists should interpret the staining results in conjunction with the patient's relevant clinical history.

Description: Prostate-specific antigen (PSA) is a single-chain glycoprotein of 237 amino acids containing approximately 8% carbohydrate. It is a serine protease produced almost exclusively by prostatic epithelial cells. Immunohistochemically PSA is expressed in the highly specialized apically-superficial layer of female and male secretory cells of the prostate gland and is readily demonstrated in adenocarcinomas of the prostate in about 99% of the cases. There is a correlation between malignancy grade and intensity of staining, high-grade carcinomas displaying a weaker expression. About 1% of poorly differentiated carcinomas have been negative for PSA. Due to the high specificity of PSA for prostatic glandular epithelium, it is very useful in identifying prostatic carcinoma in the prostate and the adjacent organs often affected by epithelial malignancies, i.e., rectum and urinary bladder. PSA may be used in a panel with NKX3.1 and Prostein, which are at least as sensitive and slightly more specific than PSA.

Supplied As: Purified antibody in Tris-HCl pH 7.4 buffer containing stabilizing proteins (including < 1% Bovine Serum Albumin) and < 0.1% ProClin by volume.

References:

1. Polascik TJ, et al. J Urol. 1999; 162:293.
2. Alanen KA, et al. Path Res Pract. 1996; 192:233-237.
3. Varma M, et al. Histopathology. 2005; 47:1-16

REF Z2710RL-A/ Z2710RS-A/ Z2710RT-A/ Z2710RP-A (1.0ml Concentrate/ 0.5ml Con./ 0.1ml Con./ 7ml Pre-dil)